

# ALP® PUMP AND SLEEVE SYSTEM

- Full range of breathable compression sleeves
- Foot, calf, full leg & pediatric SKUs
- Small and lightweight pumps (only 2.8 lbs)
- Anti-kink™ tubing to ensure air delivery



**vizient** Contract # MS8839

Currie Medical's ALP® DVT Prevention System is a non-invasive external intermittent pneumatic compression system used as a prophylaxis for reducing the incidence of Deep Vein Thrombosis (DVT). The ALP® consists of a pneumatic compression pump and an extensive range of breathable garments, including foot, calf-length, and thigh-length. The Currie ALP® System is designed with a dual focus on ease of use and patient comfort, both of which are critical components to promoting patient compliance and positive clinical outcomes.

## KEY BENEFITS

- Soft and comfortable garments increase patient compliance
- Quiet, easy to use pumps with fewer false alarms
- Universal garments fit right or left extremities or feet
- No preventative maintenance required on pumps

Currie's Uninterrupted, Continuous Sequential Fill™ technology is clinically proven to reduce the incidence of VTE.

PRODUCT ID	PRODUCT DESCRIPTION	PER CASE
ALP® 1 A	Regular CALF garment (up to 18" calf circumference)	10 pairs
ALP® 1	Regular CALF garment (up to 18" calf circumference)	10 pairs
ALP® 1 SM	Small CALF garment (up to 14" calf circumference)	10 pairs
ALP® 1P	Regular Pediatric Calf Garment (up to 13" calf circumference)	10 pairs
ALP® 1PS	Small Pediatric Calf Garment (up to 10" calf circumference)	10 pairs
ALP® 2	Large CALF garment (up to 24" calf circumference)	10 pairs
ALP® 2XL	X-Large CALF garment (up to 30" calf circumference)	10 pairs
ALP® 3	Regular THIGH garment (up to 29" thigh circumference)	10 pairs
ALP® 3 SM	Small THIGH garment (up to 25" thigh circumference)	10 pairs
ALP® 4	Large THIGH garment (up to 36" thigh circumference)	10 pairs
PVA/ALP® 1	Regular ALP® FOOT garment Universal (up to 12 EEE)	10 pairs
PVA/ALP® 2	Large ALP® FOOT garment Universal (> 12 EEE)	10 pairs
ALP® 25AK	Tubing Connector Assembly 5' Length (Anti-Kink™)	5 pairs
ALP® 25LAKP	OR Tubing Connector Assembly 10' Length (Anti-Kink™)	5 pairs