

# FLUID WARMING SYSTEMS



- Precise temperature control 95°F (35°C) - 110°F (43°C)
- Basin sizes 3.5 - 5.7 Liters
- Includes mountable stand attachment
- FDA-registered and manufactured in the USA

Northfield's heavy-duty warmer systems provide temperature control, facilitate easy drape placement and prevent drape tears and burns. Our systems help OR teams keep patients warm while maintaining a sterile field while critical procedures are performed. Northfield Medical Manufacturing is a leading provider of infection control products for procedure rooms. Health systems and surgery centers have relied on our products for over twenty years.

## KEY BENEFITS

- Access to warm irrigation, preventing unplanned hypothermia
- Adjustable, AORN-compliant temperature ranges
- Built-in caster wheels for easy movement
- Multiple basin sizes to meet your OR team's protocols

Hypothermia increases the risk of blood loss, the need for transfusions and surgical site infections.

PRODUCT ID	PRODUCT DESCRIPTION	PRODUCT DIMENSIONS
WS3035	Fluid Warming System 3.5 Liter	19.5L x 10.5W x 5H
WS3055	Fluid Warming System 5.5 Liter	19L x 12W x 5H
WS3055	Fluid Warming System 5.7 Liter	29.5L x 10.5W x 5H

Ask us about our annual warming drape plans, including **free** warming system placement!

